

SIDE RAILS - SHORT STAY

For hospitals



3/4 rails, Half-rails

REF. LSUGH2BPLTXO6, LSUGH2BLTA2, LSUGH2BPLP1/2A



A MEDJUNIC COMPANY

The side rails for AHF's hospital or ambulatory care medical beds (short stay settings) are made in our factories in France. They are designed to be **user-friendly, hygienic and lightweight while remaining robust**. Crafted from carefully selected high-quality materials, our bed rails are not only very aesthetic, but also seamlessly integrate with our AHF-made beds.

Bed rails & Product reference	Compatibility With beds	Description & Technical data	Weight
<p>Set of 2 3/4 rails REF: LSUGH2BPLTXO6 REF: LSUGH2BLTA2</p> 	<p>LSUGH2BPLTXO6 MEDIFIT</p> <p>LSUGH2BLTA2 ACTILIT</p>	<ul style="list-style-type: none"> Description: side rails covering 3/4 of the length of the bed. Two positions possible: low (non-secure) and high (secure). No tools required to remove the side rails. Rails can be positioned for an air mattress (no tools required) (MEDIFIT only) Materials: : fibreglass-reinforced plastic and epoxy-coated steel tube. Colours: cream white Dimensions: 155 x 50 x 3.6 cm (for one side rail) Standard: EN 60601-2-52 <i>(for patients having a physical size equal to or more than 146 cm, a mass equal to or more than 40 kg and a body mass index (BMI) equal to or more than 17)</i> 	<p>8 kg (for one side rail)</p>
<p>Set of 4 half-rails REF: LSUGH2BPLP1/2A</p> 	<p>ACTILIT (Half-rails included as standard with the ACTILIT bed)</p>	<ul style="list-style-type: none"> Description: 4 half-rails covering 3/4 of the length of the bed. Two positions possible: low (non-secure) and high (secure), intuitive red unlocking handle with damped descent for each rail. Options: patient and/or caregiver controls (for more information, please bed accessories data sheet – short stay) Materials: (white) polypropylene side rails with white epoxy-coated steel mechanism. Customisable laminate insert - see AHF colour chart. Dimensions: 81 x 43 x 4 cm (for one side rail) Standard: EN 60601-2-52 <i>(for patients having a physical size equal to or more than 146 cm, a mass equal to or more than 40 kg and a body mass index (BMI) equal to or more than 17)</i> 	<p>8 kg (for one side rail)</p>

Accessories available

For 3/4-length rails only

- Foam protection pads

For half-rails only

- Integrated Patient and/or Caregiver controls

Hygiene & Disinfection

Use only authorised and standardised cleaning products and disinfectants.

- Do not mix products and comply with prescribed dosages and dilutions.
- Do not apply cleaning or disinfectant products directly to the bed or accessories; apply the product to a cloth or other cleaning item and respect the drying time.
- Recommended products: simple detergents, a 97% dilution of 9 chlorometric degree bleach, products complying with standards NF EN 1040 / NF T 72-152 / NF EN 1275 / NF T 72-202, detergents containing active ingredients such as N-decyl-N, N-dimethyl-1-decanaminium chloride