

# SIDE RAILS - LONG STAY

For residential facilities for the elderly




Half-length rails, 3/4-length rails, full-length rails

REF. LSUGH2BPLTXO1/2G2, LSUGH2BPLTXO6, LSUGH2BPLPLF



A MEDJUNIC COMPANY

The side rails for AHF's hotel-style medical beds for long stay settings are made in our factories in France. They are designed to be **user-friendly and lightweight while remaining robust**. Crafted from carefully selected high-quality materials, our bed rails are not only very aesthetic, but also seamlessly integrate with our AHF-made beds.

Bed rails & Product reference	Compatibility With beds	Description & Technical data	Weight
<p>Set of 4 half-rails REF : <b>BAR1-D-G-B11 (natural wood),</b> <b>BAR1-D-G-B16 (anthracite grey)</b></p> 	<p>FLOORE VISIONA (without the connected bed rails option)</p>	<ul style="list-style-type: none"> <li>Description: half-rails covering the full length of the bed. Three positions possible: low (non-secure), intermediate and high (secure)</li> <li>Materials: grey epoxy-coated steel, aluminium and varnished solid wood. Colours (wood part): natural wood. <b>Will be available in anthracite grey in January 2025</b></li> <li>Dimensions: 92 x 53 x 7 cm (for one side rail)</li> <li>Standard: EN 60601-2-52 <i>(for patients having a physical size equal to or more than 146 cm, a mass equal to or more than 40 kg and a body mass index (BMI) equal to or more than 17)</i></li> </ul>	<p>7.2 kg (for one side rail)</p>
<p>Set of 2 3/4 rails REF: <b>LSUGH2BPLTXO6</b></p> 	<p>FLOORE VISIONA PANTOGRAPHE MEDIFEEL (without the connected bed rails option)</p>	<ul style="list-style-type: none"> <li>Description: side rails covering 3/4 of the length of the bed. Two positions possible: low (non-secure) and high (secure). No tools required to remove the side rails. Rails can be positioned for an air mattress (no tools required)</li> <li>Materials: : fibreglass-reinforced plastic and epoxy-coated steel tube. Colours: cream white. <b>Will be available in anthracite grey in January 2025</b></li> <li>Dimensions: 155 x 50 x 3.6 cm (for one side rail)</li> <li>Standard: EN 60601-2-52 <i>(for patients having a physical size equal to or more than 146 cm, a mass equal to or more than 40 kg and a body mass index (BMI) equal to or more than 17)</i></li> </ul>	<p>8 kg (for one side rail)</p>
<p>Set of 2 full-length rails REF: <b>LSUGH2BPLPLF</b></p> 	<p>FLOORE VISIONA</p>	<ul style="list-style-type: none"> <li>Description: side rail extending the full length of the bed. Four positions possible: low (non-secure), intermediate, high (secure) and very high (for air mattress)</li> <li>Material: anodised aluminium High-quality adhesive with imitation laminate texture</li> <li>Dimensions: 196 x 9.6 x 3 cm (for one side rail)</li> <li>Standard: EN 60601-2-52 <i>(for patients having a physical size equal to or more than 146 cm, a mass equal to or more than 40 kg and a body mass index (BMI) equal to or more than 17)</i></li> </ul>	<p>2.8 kg (for one side rail)</p>

## Accessories available

All bed rails

- Foam protection pad

### 3/4-rails only

- Rail extension
- Wood trim for top bar
- Laminate side panel trim

## Hygiene & Disinfection

Use only authorised and standardised cleaning products and disinfectants.

- Do not mix products and comply with prescribed dosages and dilutions.
- Do not apply cleaning or disinfectant products directly to the bed or accessories; apply the product to a cloth or other cleaning item and respect the drying time.
- Recommended products: simple detergents, a 97% dilution of 9 chlorometric degree bleach, products complying with standards NF EN 1040 / NF T 72-152 / NF EN 1275 / NF T 72-202, detergents containing active ingredients such as N-decyl-N, N-dimethyl-1-decanaminium chloride