

# BED ACCESSORIES - LONG STAY

For residential facilities for the elderly

Full list of our accessories



AHF offers a comprehensive range of accessories for its different medical bed ranges. Below is a complete list below with the corresponding bed ranges.



	Accessories and references	Compatibility with products	Description & Technical data
1	Night light <b>LSUGHVEIL</b> <b>LSUGHVEILT</b>	<b>LSUGHVEIL</b> FLOORE, VISIONA beds <b>LSUGHVEILT</b> PANTO, MEDIFEEL beds	<ul style="list-style-type: none"> <li>Description: anti-glare LED night light</li> <li>Standard: 60601-1 certificate Direct current 24V DC (maintenance-free)</li> </ul>
2	Set of 2 Foam protection pads <b>LSUGHMPBF2</b>	3/4-length rails	<ul style="list-style-type: none"> <li>Description: wraparound protection pads specially designed for AHF 3/4 side rails (for a complete bed). Easy to fit. Non-slip coating. Sold in pairs (for one bed)</li> <li>Materials: HR31030 foam, 40 mm thick Non-slip vinyl coating (grey)</li> <li>Dimensions: L 1580 x W 380 mm</li> </ul>
3	Set of 2 extensions for 3/4-rails <b>20113037</b> <b>LSUGHCBP</b>	<b>20113037</b> FLOORE, VISIONA beds <b>LSUGHCBP</b> PANTO, MEDIFEEL beds	<ul style="list-style-type: none"> <li>Description: maximise safety with these extension sections for 3/4-rails. Fasten directly onto to the bed frame at the foot of the bed. Compatible with the bed extension system. Sold in pairs (for one bed)</li> </ul>
4	Set of 2 wood trims for top bars <b>LSUGH2BPLMC</b>	3/4-length rails (long stay)	<ul style="list-style-type: none"> <li>Description: aesthetic trim fastened to the top bar of 3/4-rails, for a non-medical bedroom look. This accessory is fitted onto the side rails in the factory (option when ordering). Sold in pairs (for one bed)</li> <li>Materials: melamine</li> </ul>
5	Set of 2 laminated side panel trims <b>TOSULPAF2</b>	3/4 side rails <b>TOSULPAF2</b> FLOORE, VISIONA beds <b>TOSULPAP</b> PANTO, MEDIFEEL beds	<ul style="list-style-type: none"> <li>Description: long side panel trims positioned on the 3/4-rail to hide the bed base and give the bed a hotel-style look. Sold in pairs (for one bed)</li> <li>Materials: melamine</li> </ul>
6	Set of 4 foam protection pads <b>LSUGHMPBF1/2B</b>	Half-length rails (long stay)	<ul style="list-style-type: none"> <li>Description: wraparound protection pads specially designed for AHF long stay half-length side rails. Pads easy to fit regardless of rail position. Non-slip coating. Sold in sets of 4 (for one bed)</li> <li>Materials: HR31030 foam, 40 mm thick Non-slip vinyl coating (grey)</li> <li>Dimensions: L 970 mm x W 380 mm</li> </ul>
7	Set of 4 foam protection pads <b>LSUGHMPBFPL</b>	Full-length rails (long stay)	<ul style="list-style-type: none"> <li>Description: wraparound protection pads specially designed AHF full-length side rails. Easy to fit. Non-slip coating. Can be fitted while the bed rails are in use (4 separate modules). Sold in sets of 4 (for one bed)</li> <li>Materials: CMHR3532E foam, 15 mm thick. 6 mm spiral zipper. Vinyl non-slip coating (grey)</li> <li>Dimensions : L 1935 mm x W 300 mm</li> </ul>
8	Bedhead reading light <b>LSUGHLUMF</b>	FLOORE, VISIONA beds	<ul style="list-style-type: none"> <li>Description: map operated from the switch or the patient's control pendant. Lamp equipped with a 5.2W low-consumption LED bulb Colour temperature 3000 Kelvin. Anti-rotation system for safe use. Weight: 1.3 kg</li> <li>Dimensions: holder length 75 + 65 cm</li> </ul>
9	Compact laminate flat bed base <b>LSUGHDPICF</b> <b>LSUGHDPICV</b>	<b>LSUGHDPICF</b> FLOORE bed PANTO, MEDIFEEL beds <b>LSUGHDPICV</b> VISIONA bed	<ul style="list-style-type: none"> <li>Description: rotproof compact laminate bed base in 3 modules. Very tough, easy-care material</li> <li>Fire resistance: fire class M3</li> <li>Thickness: 6 mm, Total weight (3 modules): 10.5 kg</li> </ul>
10	Caregiver Front Panel <b>LSUGHBOITP</b>	FLOORE bed	<ul style="list-style-type: none"> <li>Description: control panel embedded at the foot of the bed for easy access by caregivers. All bed functions are available with selective locking if necessary. Weight 260 g</li> <li>Materials: Plastic casing, light grey RAL 7035 Protection IPX4</li> <li>Dimensions: 15 x 14.5 cm</li> </ul>
11	Bluetooth hand control <b>LSUGHCBF</b> <b>LSUGHCBT</b>	<b>LSUGHCBF</b> FLOORE bed (AHF Control app.) <b>LSUGHCBT</b> PANTO, MEDIFEEL beds (AHF Remote app.)	<ul style="list-style-type: none"> <li>Description: patient hand control pendant to operate all the bed functions: backrest, footrest, variable height, Trendelenburg/Anti-Trendelenburg. Bluetooth connectivity available through the AHF Control and AHF Remote apps (for IOS and Android). Selective function locking possible with the security key to be used on the back of the pendant. Protection IPX6. Weight 300 g</li> <li>Dimensions: 17 x 6 x 1 cm</li> </ul>
12	Caregiver control panel <b>LSUGHSCSTM</b>	PANTO, MEDIFEEL beds	<ul style="list-style-type: none"> <li>Description: control panel for healthcare personnel, can be hung on the footboard or 3/4-length side rail. All bed functions are available with selective locking, position memory function and emergency position.</li> <li>Materials: plastic casing, light grey RAL 7035. IPS4 protection</li> <li>Dimensions: 16 x 23 x 3 cm</li> </ul>

# BED ACCESSORIES - LONG STAY

For residential facilities for the elderly  
Full list of our accessories

# AHF

A MEDJUNIC COMPANY



Accessories and references	Compatibility with products	Description & Technical data
13 New'ee sound kit <b>LSUGHKITSONF</b>	FLOORE, VISIONA beds	<ul style="list-style-type: none"> <li>Description: sound kit for immersive sound diffusion (innovative new sound technology) directly from the headboard panel. Significantly reduces noise nuisance and improves listening quality. Consists of a transmitter box positioned near the audio source. The audio source can be a TV with a jack or RCA socket, or it can be a wireless source, such as Bluetooth, with a smartphone, for example.</li> </ul>
14 Battery <b>LSUGHBATF</b> <b>LSUGHBATAF</b>	<b>LSUGHBATF</b> FLOORE, VISIONA beds <b>LSUGHBATAF</b> PANTO, MEDIFEEL beds	<ul style="list-style-type: none"> <li>Description: optional battery useful when transporting the patient on the bed</li> <li>1.2 Ah lead gel battery IPX6 protection. Weight 1.6 kg.</li> <li>Colour RAL 7035.</li> <li>Dimensions: 11 x 12 x 8 cm</li> </ul>
15 Pillow net <b>LSUGHFILETF</b>	FLOORE, VISIONA beds PANTO, MEDIFEEL beds	<ul style="list-style-type: none"> <li>Description: net fastened under the back rest to stop the pillow from falling to the floor. Reduces musculoskeletal disorders for the nursing staff. Black polypropylene elastic net with elastic cord edging</li> <li>Dimensions: L 600 mm x W 300 mm</li> </ul>
16 Flexible pendant holder <b>LSUGHFL</b>	FLOORE, VISIONA beds PANTO, MEDIFEEL beds	<ul style="list-style-type: none"> <li>Description: holder for patient hand control pendant.</li> <li>Mounted on a flexible arm that easily fits into the anchor points provided in the bed frame</li> </ul>
17 Trapeze bar with retractable strap winder <b>LSUGHPL</b>	FLOORE, VISIONA beds PANTO, MEDIFEEL beds	<ul style="list-style-type: none"> <li>Description: trapeze bar that can be fitted onto the bed at the head end. Steel tube with strap + gripping handle. Strap retractable to adjust it to the right length between the handle and the patient. Before using, remove the trapeze handle from the storage hook.</li> <li>Trapeze bar: compliant with standard EN 60601-2-52.</li> <li>Steel tube diam. 32 mm, overall 139 x 62 cm</li> <li>Retractable strap winder: Plastic + fabric</li> </ul>
18 4-hook drip holder <b>LSUGHSTALNM</b>	FLOORE, VISIONA beds PANTO, MEDIFEEL beds	<ul style="list-style-type: none"> <li>Description: 4-hook drip holder - 2 kg max. per hook. Tube Ø25 (lower part) and Ø18 (upper part). The drip pole should not be used to move a bed, as this would risk damaging it.</li> <li>Dimensions: Closed approx. 920 mm / Open approx. 1540 mm</li> </ul>
19 Grip handle <b>20113030 B16</b>	FLOORE, VISIONA beds PANTO, MEDIFEEL beds	<ul style="list-style-type: none"> <li>Description: lifting aid with integrated mounting support. Positioned directly on the bed frame for easy grip and assistance with mobility.</li> </ul>
20 Document holder <b>LSUGHPDOC</b>	PANTO, MEDIFEEL beds	<ul style="list-style-type: none"> <li>Description: document holder (A4 size) easy to fit on the footboard panels. The edges are folded round at the bottom to stop the documents slipping out. The front is cut away to make it easier to take documents out. Weight: 2.6 kg</li> <li>Materials: steel (thickness 1.5 mm) and grey PVC coating RAL 7045</li> <li>Dimensions: L 281 mm x W 333 mm x D 129 mm</li> </ul>
21 Mattress pump bracket <b>30113060</b>	FLOORE, VISIONA beds	<ul style="list-style-type: none"> <li>Description: allows a mattress pump to be positioned on a FLOORE or VISIONA bed so that it is stable and practical for the caregiver to use. It fits onto the footboard handrail. Remark: unnecessary for the Pantographe bed, as the pump can be attached directly to the footboard with a handrail.</li> <li>Materials: galvanised steel with plastic coating</li> <li>Dimensions: wire Ø5 mm – overall size 36 x 11 cm</li> </ul>
22 Urine bag holder <b>20113011</b>	FLOORE, VISIONA beds PANTO, MEDIFEEL beds	<ul style="list-style-type: none"> <li>Description: zinc-plated steel urine bag holder. Max. load: 5 kg Weight 0.2 kg</li> <li>Materials: galvanised steel</li> <li>Dimensions: wire Ø5 mm – overall 13.5 x 5.5 cm</li> </ul>
23 Standard strap attachment <b>LSUGH1ASF</b> Large loop strap attachment <b>LSUGH1ASSPECF</b>	FLOORE, VISIONA beds PANTO, MEDIFEEL beds	<ul style="list-style-type: none"> <li>Description: from 1 to 6 strap attachments positionable on either side of the bed as needed, for a strong and effective restraint system. Resistance: 100 kg depending on force applied. Weight of one attachment: Weight 0.3 kg</li> <li>Materials: epoxy-coated steel</li> <li>Dimensions loop opening: Standard strap attachment: 7 x 1.5 cm Large loop strap attachment: 20 x 5 cm</li> </ul>
24 Set of 2 transport trolleys (for one bed) <b>LSUGHCHTRFL2</b> <b>LSUGHCHTRPAN</b>	<b>LSUGHCHTRFL2</b> FLOORE, VISIONA beds <b>LSUGHCHTRPAN</b> PANTO, MEDIFEEL beds	<ul style="list-style-type: none"> <li>Description: transport trolleys allowing the bed to be moved in a vertical position. 4 Ø100 mm castors, including 2 braked. The trolleys are easily fastened without tools to the head and foot of the bed (in the anchor points for the FLOORE bed and the bed frame for the PANTOGAPHE bed)</li> <li>Materials: epoxy-coated steel</li> <li>Dimensions: 1070 x 630 cm</li> </ul>
25 Set of 2 long wall bumpers <b>LSUGHBFF</b>	FLOORE, VISIONA beds	<ul style="list-style-type: none"> <li>Description: articulated wall bumpers 19 cm long</li> <li>Materials: aluminium, stainless steel and natural rubber NR</li> </ul>

## Hygiene & Disinfection

Use only authorised and standardised cleaning products and disinfectants.

• Do not mix products and comply with prescribed dosages and dilutions.

Do not apply cleaning or disinfectant products directly to the bed or accessories; apply the product to a cloth or other cleaning item and respect the drying time.

• Recommended products: simple detergents, a 97% dilution of 9 chlorometric degree bleach, products complying with standards NF EN 1040 / NF T 72-152 / NF EN 1275 / NF T 72-202, detergents containing active ingredients such as N-decyl-N, N-dimethyl-1-decanaminium chloride