

# TML6 Mixed, Laminate

Overbed table / U-shaped base, mixed split laminate tray

REF. GHTML6DF-S (supplied assembled, also available in kit form: see options)



The **TML6** overbed tables have been specially designed for nursing homes. Robust, ergonomic, easy to handle and space-saving, they are one of the leading products on the market. The **GHTML6DF-S** version consists of a **mechanically strong** "mixed" split laminate tray (one part fixed, one part tilting) **available in a wide choice of colours** with a **U-shaped base for a wide front opening**. AHF manufactures and offers its customers a wide range of tray types and materials (laminate, polymer and compact laminate), as well as different customisations. Our design department is also on hand to create bespoke versions to your specifications.

## Technical specifications

### TRAY

- **Mixed – split tray:**  
1 fixed part  
+ 1 tilting part

Its compact overall size is particularly suited to a 90 cm wide bed

- **Laminate**
  - > Attractive trays available in a wide choice of colours
  - > Mechanically strong
  - > Heat resistant
  - > Entirely smooth surface
  - > Rounded corners
  - > White ABS edging
- **Tray tilt**  
5-position rack mechanism

### GUARD RAIL & RETAINER

- > Guard rail fitted on the right-hand side, prevents objects falling off the tray
- > Retainer on the tilting section to prevent books or magazines slipping off

### STRUCTURE

- Structure with white epoxy-coated base + column  
Extremely sturdy and stable

### ERGONOMIC HANDLE

- Height adjustment handle located under the fixed section

### COLUMN


- The adjustable column allows use in bed or in an armchair, and adapts to the resident or patient's morphology.
- Gas spring assisted height adjustment
- With anti-pinch feature for optimum safety when raising the tray
- Aluminium column cover for a more aesthetic finish

### BASE

- Wide **U-shaped** base for excellent stability
- Base offset from the tray avoids knocking and damaging walls and the tray when moving the overbed table
- Recommended for use in an armchair
- One of the lowest bases on the market, easily fits under medical beds

### CASTORS

- 4 swivel castors (360°) dual wheel, diam. 50 mm for easy movement



## Technical data

Variable height	<b>75 cm – 106.5 cm</b>
Dimensions tilting section	<b>70 x 41 cm</b>
Dimensions fixed section	<b>20 x 41 cm</b>
Overall tray dimensions	<b>90 x 41 cm</b>
Tray thickness	<b>19 mm</b>
Base height	<b>9.5 cm</b>
Base length	<b>79 cm</b>
Base width	<b>47.5 cm</b>
Overall length	<b>90 cm</b>
Overall width	<b>47.5 cm</b>
Tilt angle	<b>65° maximum</b>
Weight when empty	<b>18.5 kg</b>
Maximum load	<b>15 kg at end of tray</b> <b>30 kg evenly distributed</b>
Materials	<b>Laminate (tray)</b> <b>Steel (structure)</b> <b>Aluminium (column cover)</b>
Dual-wheel castors	<b>360° / Ø 50 mm</b>
Fire rating (NF P 92-503)	<b>M1, non-flammable</b>

## Customisation and Accessories

See the **Colour chart** page for all the laminate colours\* available  
**More than 15 colours**

\* The laminate colours available may change (new colours, changes to shades)



### Accessories & Options available

- 2 braked castors
  - Tray edging in the same colour as the tray
  - Overbed table supplied in kit form
- REF. GHTML6DFK-S

### Hygiene & Disinfection

- Smooth surfaces, easily accessible, compatible with approved detergents. It is recommended to start cleaning at the top and work down
- Do not use abrasive substances or strippers, polish, wax, furniture cleaners or bleaching agents
- Cleaning products must be sprayed on or applied with a disposable cloth