## **SEATING ACCESSORIES**

For residential facilities for the elderly For hospitals and healthcare facilities **Full list of our accessories** 



AHF offers a comprehensive range of accessories for its armchairs Below is a complete list below with the corresponding seating ranges.



	Accessories and references	Compatibility with products	Description & Technical data
1	2 armchair castors FMFASU2G	MOMENTA and MONTARCHER armchairs	Description: 2 castors to facilitate manoeuvrability of the bedroom armchairs     Dimensions: dual wheel castor, diameter 50 mm
2	Headrest GERIART	MOMENTA and MONTARCHER armchairs	Description: headrest easily positioned on the chair back.  Non-slip on the reverse     Dimensions: 20 x 20 cm     Weight: 250 g
3	Sound kit for armchairs LSUGHKITSONFAUT	MOMENTA and MONTARCHER armchairs	Description: immersive sound diffusion directly from the armchair backrest. Reduces noise nuisance and optimises listening quality. Can be used for music therapy. Operates via a Bluetooth connection between the box fitted in the chair back and the transmitter to be connected in the bedroom

## **Hygiene & Disinfection**

- The fabrics comply with standard NF EN 15973: clean with a pH neutral soap and a soft brush. Clean stains such as pen ink, coffee, etc. immediately to avoid the stain migrating through the coating on the fabric and becoming permanent.
- Cleaning products must be sprayed on or applied with a disposable cloth
   It is also recommended to start cleaning at the top and work down

