SEATING ACCESSORIES

For residential facilities for the elderly For hospitals and healthcare facilities Full list of our accessories



AHF offers a comprehensive range of accessories for its armchairs. Below is a complete list below with the corresponding seating ranges.



	Accessories and references	Compatibility with products	Description & Technical data
1	2 armchair castors FMFASU2G	MOMENTA and MONTARCHER armchairs	 Description: 2 castors to facilitate manoeuvrability of the bedroom armchairs Dimensions: dual wheel castor, diameter 50 mm
2	Headrest GERIART	MOMENTA and MONTARCHER armchairs	 Description: headrest easily positioned on the chair back. Non-slip on the reverse Dimensions: 20 x 20 cm Weight: 250 g
3	Sound kit for armchairs LSUGHKITSONFAUT	MOMENTA and MONTARCHER armchairs	• Description: immersive sound diffusion directly from the armchair backrest. Reduces noise nuisance and optimises listening quality. Can be used for music therapy. Operates via a Bluetooth connection between the box fitted in the chair back and the transmitter to be connected in the bedroom

Hygiene & Disinfection

The fabrics comply with standard NF EN 15973: clean with a pH neutral soap and a soft brush. Clean stains such as pen ink, coffee, etc. immediately to avoid the stain migrating through the coating on the fabric and becoming permanent.
Cleaning products must be sprayed on or applied with a disposable cloth
It is also recommended to start cleaning at the top and work down

67